




Thermal Imaging

Body Temperature Measurement



With the rise in infections around the world from the Coronavirus pandemic, organisations are in need of solutions to help manage and mitigate risks.

Mitie offer solutions that can help protect your employees, visitors and ultimately, your revenue.

We partner with some of the worlds leading technology companies, meaning we can provide systems from multiple vendors. We offer a range of solutions that can be tailored to suit your needs and unique risk profile.

Our solutions range from portable, temporary deployments to fully installed permanent systems. These could be in the form of hand-held cameras your staff can use on an individual basis, or a system covering multiple entrances from one centralised location. Typical applications include health and social care, education, food retail, transportation, financial services, utilities, and other security check entrances for temporary control.



Hand Held System

A Hand-held solution that your staff can use to test individuals at the entrances to your building or randomly around your premises.



Basic System

A Basic system that can be deployed using a tripod and laptop and operated by staff at entrances. Upon detecting an individual with a high temperature, they can be diverted.



High Accuracy System


A High accuracy system that includes a temperature calibration unit. The system is also suitable for use with multiple cameras situated around site entrances and exits with monitoring from a central location.



The basis of all of these systems is thermal imaging. This technology has been used in the security industry for years and is a proven perimeter deployment that helps to secure critical sites with no loss of vision through fog or during the night.

The same technology is used to detect human body temperature. Using a single camera the accuracy is within ± 0.5 degrees centigrade. With a temperature calibration device in the field of view of the thermal camera, this accuracy is increased to ± 0.3 degrees centigrade. The temperature of someone suspected of having a fever is above 37.3 degrees. Typically the camera is calibrated to raise the alarm when it detects a temperature of 30 degrees but the calibration can lower or raise this if required.

Optimal Detection Range



Mitie have vast experience in the rapid deployment of systems and we can respond to your needs with whichever solution you require.



To speak to a representative about how we can help, contact the team on **0333 207 6606 - Option 5.**